


# Rev360 Hysteroscopy Platform

## Rev360 Hysteroscopy Platform

- » Fully rotatable stopcocks allow for easy tube orientation
- » Stryker's Speed-Lock connections allow quick assembly of hardware
- » Our 0°, 12°, 30° hysteroscopes offer unrivaled brightness, resolution, and clarity



## Rev360 IDEAL EYES Autoclavable Hysteroscopes

Part Number	Outer Diameter	Direction of View
502-729-000	2.9 mm	0°
502-729-012	2.9 mm	12°
502-729-030	2.9 mm	30°
502-740-000	4.0 mm	0°
502-740-012	4.0 mm	12°
502-740-030	4.0 mm	30°



## Rev360 Hysteroscope Sheaths

Part Number	Outer Diameter	Sheath	Use	Working Channel
502-729-051	2.9 mm	Inner	Diagnostic	N/A
502-729-052	2.9 mm	Outer	Diagnostic	N/A
502-729-061	2.9 mm	Inner	Operative	5Fr
502-729-062	2.9 mm	Outer	Operative	N/A
502-740-051	4.0 mm	Inner	Diagnostic	N/A
502-740-052	4.0 mm	Outer	Diagnostic	N/A
502-740-061	4.0 mm	Inner	Operative	7Fr
502-740-062	4.0 mm	Outer	Operative	N/A

## Rev360 Hysteroscope Obturators

Part Number	Sheath Diameter	Sheath Description
502-729-060	2.9 mm	Operative
502-729-050	2.9 mm	Diagnostic
502-740-060	4.0mm	Operative
502-740-050	4.0 mm	Diagnostic

\*This does not replace the operative technique, please refer to the operative technique for all procedures. \*\*Not drawn to scale.

# Rev360 Instrumentation



## Rev360 5Fr Semi-Flexible Instruments

For Use with 2.9 mm or 4.0 mm  
Operative Sheath

- 502-555-100 Semi-Flexible (Single Action) Grasper, 5Fr, 41 cm
- 502-555-200 Semi-Flexible Biopsy Punch, 5Fr, 41 cm
- 502-555-300 Semi-Flexible Scissors, 5Fr, 41 cm
- 502-555-400 Semi-Flexible Biopsy Cup, 5Fr, 41 cm
- 502-555-500 Semi-Flexible (Double Action) Grasper, 5Fr, 41 cm

## Rev360 7Fr Semi-Flexible Instruments

For Use with 4.0 mm Operative Sheath

- 502-777-100 Semi-Flexible (Single Action) Grasper, 7Fr, 41 cm
- 502-777-200 Semi-Flexible Biopsy Punch, 7Fr, 41 cm
- 502-777-300 Semi-Flexible Scissors, 7Fr, 41 cm
- 502-777-400 Semi-Flexible Biopsy Cup, 7Fr, 41 cm
- 502-777-500 Semi-Flexible (Double Action) Grasper, 7Fr, 41 cm

## Rev360 HTA<sup>®</sup> Adapter (Hydro ThermAblator<sup>®</sup>)


- 502-023-029 HTA Adapter for 2.9 mm Hysteroscope

## Rev360 Accessories


- 502-880-205 Port Seal (10/box)
- 502-880-650 Reusable Stopcock

\*This does not replace the operative technique, please refer to the operative technique for all procedures. \*\*Not drawn to scale.

# Rev360 22/24Fr Resection Platform



## REV360 22/24FR RESECTION PLATFORM

<b>Rev360 Resectoscope</b>	<b>502-729-012</b> 2.9 mm IDEAL EYES Autoclavable Hysteroscope
<b>Rev360 Resectoscope Sheaths</b>	<b>502-990-322</b> 22Fr Resectoscope Standard Flow Sheath (Inner) <b>502-990-424</b> 22/24Fr Continuous Flow Resectoscope Sheath Set (Inner/Outer)
<b>Rev360 Resectoscope Obturator</b>	<b>502-990-414</b> 22/24Fr Resectoscope Standard Obturator
<b>Rev360 Resectoscope Working Element</b>	<b>502-990-401</b> 22/24Fr Working Element—Passive
<b>Rev360 Resectoscope Electrodes</b> (Sterile/Disposable) For use with working element part #502-990-401 only (Quantity: box of 6)	<b>504-990-314</b> Rev360 22Fr 30° Cutting Loop Electrode, .014 <b>504-990-200</b> Rev360 22Fr Roller Ball Electrode <b>504-990-600</b> Rev360 22Fr Dimpled Vaporizing Roller Ball Electrode <b>504-990-100</b> Rev360 22Fr Gold Coagulating Electrode <b>504-990-500</b> Rev360 22Fr 0° (Straight) Cutting Loop Electrode, .014
	<b>Sterilization Trays</b> <b>233-032-880</b> Urology/Gynecology Accessory Tray (Single Tier) <b>233-032-881</b> Urology/Gynecology Tray (Two Tiers)

\*This does not replace the operative technique, please refer to the operative technique for all procedures. \*\*Not drawn to scale.

# 24/26FR 4 MM Resection PProducts



**Resectoscope**



**Resectoscope Sheaths**



**Resectoscope Working Elements**

<b>Resectoscope</b>	<b>502-990-030</b> 4.0 mm Scope use with 24/26Fr Resectoscope
<b>Resectoscope Sheaths</b>	<b>502-880-324</b> 24Fr Standard Flow Resectoscope Sheath (Inner) <b>502-880-426</b> 24/26Fr Continuous Flow Resectoscope Sheath Set (Inner/Outer)
<b>Resectoscope Working Elements</b>	<b>502-880-401</b> Working Element—Passive (Iglesias) <b>502-880-402</b> Working Element—Active (Baumrucker)

<b>Resectoscope Obturators</b>	<b>502-880-001</b> 24/26Fr Resectoscope Standard Obturator <b>502-880-002</b> 24/26Fr Resectoscope Visual Obturator
--------------------------------	--

<p><b>Resectoscope Electrodes</b> (Sterile/Disposable) For use with 24/26Fr Resectoscope, Working Element part #502-880-401 and #502-880-402 (Quantity: box of 6)</p>	<p><b>504-880-414</b> 24Fr 30° Cutting Loop Electrode, .014 <b>504-880-314</b> 24Fr 90° Cutting Loop Electrode, .014 <b>504-880-200</b> 24Fr Roller Ball Electrode <b>504-880-600</b> Dimpled Vaporizing Roller Ball Electrode <b>504-880-100</b> Gold Coagulating Electrode <b>504-880-500</b> 0° (Straight) Cutting Loop Electrode <b>504-880-800</b> Collings Knife Electrode</p>
---	--

\*This does not replace the operative technique, please refer to the operative technique for all procedures. \*\*Not drawn to scale.

# Resection Accessories

## Resectoscope Cables

- 502-990-100** Aesculap/BERCHTOLD Cable
- 502-990-200** Storz/Erbe Cable
- 502-990-300** Universal/Valleylab Cable

## Rigid Optical Instruments

For use with 24/26Fr Resectoscope only



- 502-880-210** Rigid Optical Biopsy Forcep 0°, 12°
- 502-880-211** Rigid Optical Grasping Forcep 0°, 12°
- 502-880-213** Rigid Optical Stone Crushing Forcep 30°
- 502-880-214** Rigid Optical Biopsy Forcep, 70°
- 502-880-215** Rigid Optical Biopsy Forcep, 30°
- 502-880-230** Rigid Optical Stent Removal Forcep 0°, 12°

## Resectoscope Bridges and Adapters



- 502-880-006** Bladder Syringe Adapter, Inner Sheath
- 502-880-007** Bladder Syringe Adapter, Outer Sheath
- 502-880-600** Ellik Evacuator
- 502-880-602** Ellik Adapter, Inner Sheath
- 502-880-603** Ellik Adapter, Outer Sheath
- 502-880-005** Laser Bridge
- 502-880-208** Rigid Optical Forcep Adapter, Outer Sheath
- 502-880-209** Rigid Optical Forcep Adapter, Inner Sheath

## Resectoscope Accessories


- 502-880-205** Port Seal (10/box)
- 502-880-650** Reusable Stopcock

\*This does not replace the operative technique, please refer to the operative technique for all procedures. \*\*Not drawn to scale.

# Fluid Management

## FLUID SAFE Fluid Management System

- » Eliminates procedure interruption with automatic fluid bag recognition
- » Programmable flow and pressure rates per surgeon preference
- » Patient safety features including automatic deficit measurement



## FLUID MANAGEMENT

<b>FLUID SAFE Hysteroscopy Pump</b>	<b>502-000-000</b>	FLUID SAFE Fluid Management System (3 Piece Kit) Includes monitor, pump console, and wheel base
	<b>502-000-004</b>	Fluid Safe Pump incl. IFU   EN,NL,FR,DE, PL and Power cord EURO 2m
	<b>502-000-005</b>	Fluid Safe Pump incl. Manual EN,ES,IT,PT,GR and Power cord EURO 2m
	<b>502-000-006</b>	Fluid Safe Pump incl. Manual EN,NO,SV,DA,FI and Power cord EURO 2m
	<b>502-000-002</b>	Fluid Safe Weighing System
	<b>502-000-003</b>	Fluid Safe Roller Base
	<b>Z0176-01</b>	Power Supply Cord UK, 2 M
	<b>Z0177-01</b>	Power Supply Cord CH, 2 M
	<b>Z1409-39</b>	In-/Outflow tube set, disposable (box of 10 pcs.)
	<b>Z1420-39</b>	Tube set pump PVC, disposable (box of 10 pcs.)
	<b>Z1421-39</b>	Tube set pump PVC , safe lock, disposable (box of 10 pcs.)
	<b>Z1428-39</b>	Tube set pump SIL, reusable
	<b>Z1406-39</b>	Outflow tube set, disposable (box of 10 pcs.)
	<b>Z0176-01</b>	Power Supply Cord UK, 2 M
	<b>Z0177-01</b>	Power Supply Cord CH, 2 M
	<b>Z0554-01</b>	Container Bemis for pump 2 pcs.
	<b>Z0555-01</b>	Container Bemis for container scale 36 pcs.
	<b>Z0553-01</b>	Coil of wire for container Bemis
	<b>2200098</b>	Container adapter for VacSax (box of 5 pcs.)
	<b>2200097</b>	Container adapter for Medela
	<b>2203086</b>	Container adapter for Abbott

\*This does not replace the operative technique, please refer to the operative technique for all procedures. \*\*Not drawn to scale.